
# 106-101 Lozenge Anti-Skate Stud



#### Description

Manufactured from 316 grade steel. Satin polish finish as standard. Constructed from non corrosive material. Decorative feature anti-skating device. 2no. M8 studs to be resin fixed to concrete, stone or timber.

#### **Dimensions**


25mmø


80mm Overall length + 12.5mm height

4mm Internal radius

### **Options**

Specify required internal radius to suit edge of base product. 304 Grade stainless steel Bright polish finish Bespoke sizing Full installation service








## 106-101 Lozenge Anti-Skate Stud

### Installation



### Installation

- 1. Drill 2no. holes 10mm in diameter to the depths shown above.
- 2. Remove dust and debris caused by drilling.
- 3. Following the product guidelines heat the Polymer Modified Adhesive and pour a sufficient amount to fill the drilled hole.
- 4. Push the stud into the hole firmly ensuring an edge seal around the entire stud
- 5. Allow a few minutes for the Polymer Modified Adhesive to set.

### Cleaning & Maintenance

This product is designed for all weather exposure requiring only minimal maintenance as detailed below.

<u>Stainless Steel:</u> Manufactured from premium grade stainless steel. Maintenance free, inspect fixing every 6 months. Add additional Adhesive where required.



# 106-101 Lozenge Anti-Skate Stud

